Sinus Lift & Implantology Bone Instruments

Sinus Lift Instruments, Osteotomes, Bone Spreaders, Bone Condensers, Hammers/Mallets

en-English

Attention

The operator/user must carefully read and understand this manual thoroughly to keep the product performance durable and reliable for defined life expectancy. After opening the packaging, first of all it is necessary to check the component(s) against the standard configuration. Check that these are all present and in perfect condition.

Intended Use:

These instruments are intended to be used for various oral and maxillofacial surgical procedures, including sinus membrane elevation, bone preparation, condensation, expansion, and implant site optimization. They facilitate precise manipulation of hard tissues while minimizing trauma and preserving surrounding anatomical structures.

The Instrument Specific Intended Use along with their variants is given below:

Sinus Lift Instruments:

Sinus lift instruments are specifically designed for the atraumatic elevation of the Schneiderian membrane during maxillary sinus augmentation procedures. They assist in creating sufficient bone height for the placement of dental implants in the posterior maxilla.

Osteotomes:

These are intended to be used for the controlled condensation, expansion, and preparation of bone in the maxilla or mandible to optimize implant site density and stability. They facilitate the creation of osteotomies with minimal trauma, enhancing primary implant stability.

Bone Spreaders:

Bone spreaders are designed to aid in the separation and expansion of bone structures, particularly the alveolar ridge, to accommodate dental implants. They enable the widening of narrow ridges with minimal risk of fracture, allowing for the placement of implants in anatomically challenging sites.

Bone Condensers:

Bone condensers are intended to be used for compacting bone graft materials into osteotomy sites, ridges, or sinus windows during grafting procedures. They ensure optimal density and stability of the graft material, promoting successful osseointegration of dental implant.

Hammers/Mallets

Hammers and mallets are intended for the controlled application of force during dental and oral surgical procedures. They assist in bone manipulation, osteotomy, and the placement of other surgical instruments, ensuring precise and controlled force application.

 Our Dental Bone & Hard Tissue Surgical Instruments are intended to be used by qualified dental professionals only in a professional dental facility.

Patient population:

These instruments are intended to be used on adult patients undergoing oral and maxillofacial surgical interventions. They are primarily used in individuals requiring sinus lifts, bone grafting, implant site preparation, or ridge expansion. The appropriate instrument type and size must be selected based on the patient's oral anatomy, surgical site, and procedural requirements.

For Available Configurations/Models please visit the Implantology section in our Catalogue

Features

Sinus Lift Instruments:

- Unique curvatures and angulations for precise and atraumatic membrane elevation.
- Rounded or spoon-shaped tips to prevent tearing of the Schneiderian membrane.
- Varying lengths and diameters for better anatomical access.
- Smooth, polished surfaces reduce friction and minimize tissue damage.

Osteotomes:

- Designed for controlled bone condensation and ridge expansion.
- High-precision stainless steel or titanium construction for durability.
- Ergonomic anti-slip handles for better grip and control.

Bone Spreaders:

- Allow gradual and controlled ridge expansion.
- Made from high-grade stainless steel, resistant to corrosion.
- Designed for predictable bone spreading without excessive trauma.

Bone Condensers:

- Compact bone graft materials efficiently in sinus lift and ridge augmentation procedures.
- Feature precise tip designs for controlled bone compaction.
- Dual-ended designs for versatility in bone manipulation.

Hammers/Mallets:

- Ergonomic non-slip grips for controlled force application.
- Shock-absorbing mechanisms reduce operator fatigue.
- Designed for use in bone manipulation and osteotomy procedures.

Operating instructions

 These instruments are non-sterile and must be cleaned and sterilized before first use.

 Always check if the reprocessed instruments are in perfect condition and do not use the instruments if any kind of deterioration, damage or deformation is observed.

Sinus Lift Instruments:

- Select the appropriate sinus lift instrument for the patient's anatomy.
- Carefully introduce the instrument into the surgical site and gently elevate the sinus membrane using a controlled lifting motion.
- Ensure smooth and continuous elevation to avoid perforation.

Osteotomes

- Position the osteotome at the osteotomy site and apply controlled pressure.
- Use a mallet for gradual bone expansion as needed.
- Monitor bone density changes and adjust force application accordingly.

Bone Spreaders:

- Insert the spreader into the prepared site and apply gradual force to widen the ridge.
- Ensure that expansion is done incrementally to minimize fracture risk.

Bone Condensers:

- Use gentle, tapping motions to compact bone graft material.
- Alternate compaction and assessment for optimal graft stability.

Hammers/Mallets:

- Use in combination with osteotomes or chisels for controlled bone manipulation.
- Apply precise and controlled force to avoid excessive trauma to surrounding tissues.

 Avoid excessive mechanical impact on delicate instruments like sinus lift tools and condensers. Use hammers and mallets with controlled force to avoid excessive trauma to bone structures.

Adverse Event

Any serious incident that has occurred concerning the device should be reported to the manufacturer and to the respective competent authority of the Member State in which the incident has occurred.

Cleaning & Sterilization

Cleaning

 Clean the device with care after each application to avoid the drying of contaminants.

Initial Cleaning

1. Remove any visible debris from the instrument using moisten a cotton swab or cloth with isopropyl or ethyl alcohol and gently wiping the surface of the product.
2. Use a soft brush and enzymatic cleaner to remove contaminants thoroughly.

Manual Cleaning

3. Prepare a solution of neutral pH enzymatic detergent according to the manufacturer's instructions.
4. Immerse in a pre-soak enzymatic cleaner solution
5. Thoroughly clean all surfaces of the Instrument using a soft brush or cloth. Pay close attention to any crevices or hard-to-reach areas.
6. Rinse thoroughly under lukewarm running tap water for a minimum of 30 seconds to remove all traces of detergent.

Automated Cleaning

1. Always clean in an appropriate box or cassette and ensure that they are properly positioned within the unit.
2. Always follow the equipment manufacturer's instructions and ensure the instruments are compatible with the cleaning system.
3. Visually inspect for cleanliness, If any contamination is visible, repeat the cleaning steps until required cleanliness is achieved.
4. Dry with clean, lint free cloth or filtered compressed air until there is no visible moisture. Follow STEAM STERILIZATION PROCEDURE.

 Ultrasonic cleaning may be performed in accordance with hospital or facility protocols. Ensure proper rinsing and drying after the process to prevent residue buildup.

Sterilization

Note: Ensure that the instruments are completely dry before sterilization.

Steam Autoclaving

1. Place the instrument in an approved sterilization pouch or wrap.
2. Sterilize using a steam autoclave following a validated sterilization cycle according to ISO 17665-1 and your facility's protocols. A typical recommended cycle is 132-134°C (270-273°F) for 4 minutes.
3. Verify no residual moisture remains after sterilization to prevent corrosion and ensure longevity.
4. Allow the instruments to cool and dry completely before handling. Store in a clean, dry, and controlled environment until ready for use.
5. Follow the pouch/wrap manufacturer instructions for storage conditions and maximum storage time.

Important Notes

- Always Inspect the instruments for any signs of deterioration, breakage and bent etc. after each sterilization cycle. Discard if damage is observed.
- Always follow your facility's infection control protocols.
- Refer to your autoclave manufacturer's instructions for proper operation and maintenance. The Sterilization equipment must be validated by the hospital and or sterilization equipment manufacturers.
- The above-mentioned cleaning and sterilization guidelines, provided by manufacturers, are intended as procedures compatible with specific materials. Sterilization must be performed according to the Hospital/Clinic approved protocol.
- The responsibility for ensuring proper sterilization of instruments lies with the user. The effectiveness of sterilization depends on validated cleaning, packaging, and sterilization procedures carried out at the facility.

Disposal

Our Sinus Lift and Bone Implantology Instruments, Osteotomes, Bone Spreaders, Bone Condensers, Hammers/Mallets contain no hazardous materials. However, these must be cleaned/sterilized as per defined procedures before disposition. Please follow the local and national regulations or healthcare facilities' defined disposal and waste management policies.

Disclaimer

The products must be used, reprocessed, and maintained strictly in accordance with the instructions provided above. Any deviation from these guidelines by the dental professional or

user is undertaken at the user's sole risk. The Manufacturer will not accept any requests for refunds or exchanges under warranty for products that have not been handled and reprocessed in compliance with these instructions.

Explanation of utilized symbols

	Medical Device Compliant with Regulation (EU) 2017/745
	Consult instructions for use
	Medical Device
	Product Code/Catalogue Number
	Model/Variant
	Lot Number
	The Device is supplied non-sterile.
	Sterilizable in a steam sterilizer (autoclave) at 132-134°C (270-273°F)
	Manufacturing Date
	Manufacturer
	Authorized representative in the European Community
	Importer
	Keep in a cool, dry place
	Caution, read instructions

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