# Ecosite Elements

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English

## Product description

Ecosite Elements is a light-curing, esthetic, radiopaque nanohybrid composite system for the anterior and posterior region. The range of highly viscous composites includes 5 VITA shades (B1, A2, A3, A3.5 and A4) and 4 enamel shades (Enamel Bleach = EB, Enamel Light = EL, Enamel Medium = EM and Enamel Dark = ED) in different opacities. This, combined with their unique handling and excellent polishing properties contributes to the creation of highly-aesthetic restorations.

The shades OA2 (Opak A2) and INC (Incisal) are radiopaque composites with optimum flow composites. The shade OA2 is recommended for the veneering of slightly discolored areas and for creating ultimate tooth opacity. The shade INC is particularly suitable for an esthetic incisal edge design thanks to its high transparency. Both shades are suitable for fissure sealing.

The shades W (White) and B (Brown) are flowable composites for the customization of composite restorations. The shade W is used for the reproduction of white spots and as an opaque covering layer. The shade B is an ideal choice for fissures and the replication of brown spots.

## Intended purpose

### Ecosite Elements in the colors B1, A2, A3, A3.5, A4, EB, EL, EM, ED

- Restorations of all cavity classes of every dentition
- Direct inlays, onlays and veneers (vestibular and occlusal)
- Core build-up
- Splinting/blocking of loose teeth
- Repairs

### Ecosite Elements in the colors OA2 and INC

- Restorations of all cavity classes III-V of every dentition
- Liner
- Fissure sealing
- Splinting/blocking of loose teeth
- Repairs

### Ecosite Elements in the colors W and B:

- Customization of composite restorations

## Limitations of use

Application of the material is contraindicated if dry isolation or the recommended application technique is not possible.

## Indications

- Loss of tooth substance due to caries, fractures, tooth abrasion or development disorders
- Diastema, tooth discoloration or slight misalignment of the anterior teeth
- Dental or periodontal diseases, e.g. trauma and periodontitis, which require the splinting for the affected teeth
- Primary and secondary prevention of caries

## Contraindications

- Do not use directly on the exposed pulp.
- Do not use the material where there are known allergies to any of the components or contact allergies.

## Patient target group

People treated in the course of a dental procedure.

## Intended users

Dentist

## Notes for use

- For optimum polymerization, position the light-curing unit as close to the restorative material as possible and observe the recommendations relating to the light-curing time and the light intensity (see table). Observe the manufacturer's instructions for the light-curing unit used.
- Follow the manufacturer's instructions for other materials that are used with Ecosite Elements.
- If the material is applied into the patient's mouth with a single use applicator, the single use applicator must only be used on one patient for hygiene reasons.

## Recommended use

(See the end of the instructions for use)

1. Clean the tooth. Remove all cleaning residue with water spray. Determine tooth shade with the Ecosite Shade Guide Pure, the Ecosite Shade Guide Layer or the VITA tooth shade. Select suitable Ecosite Elements shade or shade combinations. Either choose a monochromatic Pure shade or combine it with a layer shade according to the outline. The specification corresponds to the use of a 0.5 mm-thick Ecosite Elements Layer:

**Note:** DMG recommends the use of a rubber dam.

2. Prepare the cavity in accordance with the general rules for adhesive procedures. Bevel the enamel edges. Ensure that the areas close to the pulp are protected.
3. Condition the dental enamel and apply the adhesion primer according to the manufacturer's instructions for the material used. DMG advises using the universal adhesion primer Ecosite Bond.
4. If necessary apply an underfilling of low-viscosity composite or compomer. We recommend using a matrix.
5. Insert Ecosite Elements with a maximum layer thickness of 2 mm into the cavity. Apply taller restorations one layer at a time. Light-cure each layer separately.

Light intensity	> 1,000 mW/cm <sup>2</sup>	> 600 mW/cm <sup>2</sup>
Radiation time	20 s	40 s
Max. layer thickness	2 mm	2 mm
Wavelength	420 – 480 nm	420 – 480 nm

**Note:** Additional buccal or lingual light-curing is advisable. Remove metal matrices beforehand.

- Removal of excess material, contouring and fine processing are done using diamond finishing tools and flexible polishing discs of decreasing grain size, for example. Ecosite Elements can be given a high shine using polishing paste together with polishing cups or discs.

#### Customization of composite restorations with W (White) and B (Brown)

- Select the color required for the desired color shade.
- Pour a sufficient amount of the selected shade onto a mixing pallet.
- Apply the material to the composite with a fine brush in a thickness of  $\leq 0.1$  mm. Light-cure each layer separately.

**Note:** Always cover the customization layer with a composite layer.

**Note:** Mix the customization shades with one another (to obtain new tones), but never with other products.

Light intensity	> 1,000 mW/cm <sup>2</sup>	> 600 mW/cm <sup>2</sup>
Radiation time	20 s	40 s
Max. layer thickness	0.1 mm	0.1 mm
Wavelength	420 – 480 nm	420 – 480 nm

## Disinfection of the Safetip-Dispenser Pro

The Safetip-Dispenser Pro is reusable after disinfection, although it should be replaced if visible signs of wear appear.

The Safetip-Dispenser Pro is suitable for steam sterilization in autoclaves. The procedure must be carried out at 134°C/3 bar for 5 minutes.

### Interactions

- Materials containing eugenol, moisture and oily air may inhibit polymerization at the contact area. Avoid this.
- The use of oral rinses, plaque indicators and chlorhexidine may lead to discoloration.

### Residual risks/side effects

There are no known side effects to date. The residual risk of hypersensitivity to components of the material cannot be ruled out.

### Warnings / precautions

- Keep out of the reach of children!
- Avoid contact with the skin! In the event of accidental skin contact immediately wash the affected area thoroughly with soap and water.
- Avoid eye contact! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Reuse: Standard surface disinfection as required. No specific maintenance required.
- Serious incidents involving this product must be reported to the manufacturer and to the responsible registration authorities.

### Performance characteristics of the device

Ecosite Elements in the shades B1, A2, A3, A3.5, A4, EB, EL, EM, ED complies with DIN EN ISO 4049: 2009 (Type 1, Class 2, Group 1); polymerization depth:  $\leq 2$  mm with a light-curing time of 40 seconds with a light intensity of  $\geq 600$  mW/cm<sup>2</sup> and a light-curing time of 20 seconds with a light intensity of  $\geq 1000$  mW/cm<sup>2</sup>

Ecosite Elements in shades OA2 and INC correspond to DIN EN ISO 4049:2009 (Type 2, Class 2, Group 1); polymerization depth  $\geq 2$  mm with a light-curing time of 40 seconds with a light intensity of  $\geq 600$  mW/cm<sup>2</sup> and a light-curing time of 20 seconds with a light intensity of  $\geq 1000$  mW/cm<sup>2</sup>

### Storage / disposal

- Store in a dry place (2-25°C/36-77°F) at room temperature.
- Dispose of in accordance with national regulations.

### Composition

#### Ecosite Elements B1, A2, A3, A3.5, A4, EB, EL, EM, ED:

(Barium) dental glass (approx. 81 %) in a resin matrix based on EBPADMA and Bis-GMA (approx. 19%). Pigments, additives and catalysts.

Filler content: approx. 81 wt.% = 65 vol.%. The inorganic filler particle size varies between 0.02 and 0.7  $\mu$ m.

#### Ecosite Elements OA2, INC, W and B:

(Barium) dental glass (approx. 65 %) in a resin matrix based on EBPADMA, Bis-GMA and TEDMA (about 34%)... Pigments, additives and catalysts.

Filler content: approx. 65 wt.% = 38 vol.%. The inorganic filler particle size varies between 0.02 and 0.7  $\mu$ m.

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