Dental Impression and Anesthesia Instruments

Impression Trays, Dental Syringes(Without needles), Syringe Holders

ErgoDenta ApS

Hørskættø 14, 2630 Taastrup, Denmark.

Tel: +4555241100

Email: info@ergodenta.dk

Web: www.ergodenta.dk

en-English

Attention

The operator/user must carefully read and understand this manual thoroughly to keep the product performance durable and reliable for defined life expectancy. After opening the packaging, first of all it is necessary to check the component(s) against the standard configuration. Check that these are all present and in perfect condition.

Intended Use:

These instruments are intended to be used for various dental procedures, including impression-taking, anesthetic administration, and syringe organization.

The Instrument Specific Intended Use along with their variants is given below:

Impression Trays:

Impression Trays are intended for holding and supporting impression materials to create accurate negative molds of the dental arches, teeth, and surrounding oral structures. They assist in the fabrication of dental prosthetics, restorations, and orthodontic appliances.

Variants:

Perforated: It has small holes across its surface, allowing for better retention of the impression material.

Non-Perforated:

It has a smooth surface and is typically used with materials that adhere well on their own. It is used in cases where a more precise impression is needed.

Dental Syringes (Without needle):

These are intended to be used for the controlled administration of anesthetic solutions, irrigation fluids, or other dental materials during dental procedures. They assist in precise delivery and application to support clinical treatments.

Variants:

- **Pen Model Intraligamental Syringe** – Designed for controlled intraligamental anesthesia administration.
- **Cartridge Syringe (Metric)** – Compatible with standard anesthetic cartridges for precise injection.
- **Self-Aspirating Syringe** – Equipped with an aspiration function to enhance safety by preventing intravascular injection.

Syringe Holders:

Syringe holders are intended to securely hold dental syringes or needle sheaths during procedures, minimizing the risk of accidental needle sticks and maintaining an organized workspace. They assist dental professionals by providing a stable platform for syringe placement, enhancing safety and efficiency.

 Our Impression Trays, Dental Syringes, Syringe Holders Instruments are intended to be used by qualified dental professionals only in a professional dental facility.

Patient population:

These instruments are intended to be used on adult patients undergoing restorative dentistry, prosthodontics, orthodontic treatments, oral surgery, periodontal procedures, and general dental anesthesia applications

For Available Configurations/Models please visit the Diverse section in our Catalogue

Features

Impression Trays:

- Engineered to provide a rigid and stable support structure, minimizing distortion of the impression material during setting and removal, ensuring accurate replication of dental arches and oral anatomy.
- Available in standardized sizes (small, medium, large) and configurations (full-arch, quadrant, sectional) to accommodate diverse anatomical variations and procedural requirements.
- Fabricated from high-grade, biocompatible stainless steel, ensuring sterilizability, durability, and compliance with infection control protocols.
- Designed with optimized peripheral borders to support proper extension of the impression material.

Dental Syringes:

- Designed to deliver local anesthesia with high precision, ensuring targeted numbing of specific oral tissues while minimizing patient discomfort.
- Equipped with an aspirating function, allowing dentists to check for blood return and avoid injecting anesthetic into blood vessels.
- Lightweight and ergonomically designed for easy handling, providing better control and reducing hand fatigue during procedures.
- Compatible with standard anesthetic cartridges and needles, ensuring seamless integration into dental practices.
- Made from durable, autoclavable materials, allowing for repeated use while maintaining strict hygiene standards.

Syringe Holder:

- Provides a stable and secure holding mechanism for various types of dental syringes, ensuring easy accessibility during procedures.
- Designed to accommodate different syringe models, including pen model intraligamental syringes, cartridge syringes, metric syringes, and self-aspirating syringes.
- Constructed from high-quality, autoclavable materials, ensuring sterilizability and compliance with infection control standards.
- Ergonomically designed for ease of use, reducing the risk of accidental drops or contamination.
- Compact and lightweight, allowing convenient placement within the clinical workspace without occupying excessive space.

Operating instructions

⚠ These instruments are non-sterile and must be cleaned and sterilized before first use.

⚠ Always check if the reprocessed instruments are in perfect condition and do not use the instruments if any kind of deterioration, damage or deformation is observed.

Impression Trays:

- Select the appropriate size and type of impression tray based on the patient's oral anatomy and procedure requirements.
- Load the impression material evenly into the tray, avoiding air pockets.
- Position the tray in the patient's mouth, ensuring full coverage of the required dental arch.
- Apply gentle, even pressure to seat the tray properly without causing discomfort to the patient.
- Hold the tray in position until the impression material has set, then carefully remove it in a single, smooth motion to prevent distortion.

Dental Syringes:

- Select the appropriate syringe type and ensure compatibility with the intended anesthetic or fluid.
- Load the anesthetic cartridge or fluid into the syringe according to manufacturer instructions.
- Attach a sterile needle by twisting it securely onto the syringe hub.
- Check for proper function, ensuring smooth aspiration and injection capability.
- Administer the anesthetic or fluid using the correct technique based on the procedural requirement, then safely remove the needle and dispose of it in a designated sharps container.

Syringe Holder:

- Ensure the syringe holder is clean, sterilized, and free of contaminants before use.
- Select a stable and accessible placement within the clinical workspace.
- Insert the syringe into the designated slot or holding mechanism, ensuring a secure fit.
- Remove the syringe smoothly when needed, avoiding unnecessary force.

⚠ Avoid excessive force or improper technique when using these instruments, as this can lead to patient discomfort, breakage, or procedural inefficiencies.

Adverse Event

Any serious incident that has occurred concerning the device should be reported to the manufacturer and to the respective competent authority of the Member State in which the incident has occurred.

Cleaning & Sterilization

Cleaning

Clean the device with care after each application to avoid the drying of contaminants.

Initial Cleaning

1. Remove any visible debris from the instrument using a moisten cotton swab or cloth with isopropyl or ethyl alcohol and gently wiping the surface of the product.
2. Impression trays should have residual impression material removed before it hardens using a soft brush and running water.

Manual Cleaning

3. Prepare a solution of neutral pH enzymatic detergent according to the manufacturer's instructions.
4. Immerse in a pre-soak enzymatic cleaner solution
5. Thoroughly clean all surfaces of the Instrument using a soft brush or cloth. Pay close attention to any crevices or hard-to-reach areas.
6. Rinse thoroughly under lukewarm running tap water for a minimum of 30 seconds to remove all traces of detergent.

Automated Cleaning

1. Always clean in an appropriate box or cassette and ensure that they are properly positioned within the unit.
2. Always follow the equipment manufacturer's instructions and ensure the instruments are compatible with the cleaning system.
3. Visually inspect for cleanliness, If any contamination is visible, repeat the cleaning steps until required cleanliness is achieved.

4. Dry with clean, lint free cloth or filtered compressed air until there is no visible moisture. Follow STEAM STERILIZATION PROCEDURE.

Sterilization

Note: Ensure that the instruments are completely dry before sterilization.

Steam Autoclaving

1. Place the instrument in an approved sterilization pouch or wrap.
2. Sterilize using a steam autoclave following a validated sterilization cycle according to ISO 17665-1 and your facility's protocols. A typical recommended cycle is 132-134°C (270-273°F) for 4 minutes.
3. Verify no residual moisture remains after sterilization to prevent corrosion and ensure longevity.
4. Allow the instruments to cool and dry completely before handling. Store in a clean, dry, and controlled environment until ready for use.
5. Follow the pouch/wrap manufacturer instructions for storage conditions and maximum storage time.

Important Notes

- Always Inspect the instruments for any signs of deterioration, breakage and bent etc. after each sterilization cycle. Discard if damage is observed.
- Always follow your facility's infection control protocols.
- Refer to your autoclave manufacturer's instructions for proper operation and maintenance. The Sterilization equipment must be validated by the hospital and or sterilization equipment manufacturers.
- The above-mentioned cleaning and sterilization guidelines, provided by manufacturers, are intended as procedures compatible with specific materials. Sterilization must be performed according to the Hospital/Clinic approved protocol.
- The responsibility for ensuring proper sterilization of instruments lies with the user. The effectiveness of sterilization depends on validated cleaning, packaging, and sterilization procedures carried out at the facility.

Disposal

Our Impression Trays, Dental Syringes, Syringe Holders contain no hazardous materials. However, these must be cleaned/sterilized as per defined procedures before disposition. Please follow the

local and national regulations or healthcare facilities' defined disposal and waste management policies.

Disclaimer

The products must be used, reprocessed, and maintained strictly in accordance with the instructions provided above. Any deviation from these guidelines by the dental professional or user is undertaken at the user's sole risk. The Manufacturer will not accept any requests for refunds or exchanges under warranty for products that have not been handled and reprocessed in compliance with these instructions.

Explanation of utilized symbols

	Medical Device Compliant with Regulation (EU) 2017/745
	Consult instructions for use
	Medical Device
	Product Code/Catalogue Number
	Model/Variant
	Lot Number
	The Device is supplied non-sterile.
	Sterilizable in a steam sterilizer (autoclave) at 132-134°C (270-273°F)
	Manufacturing Date
	Manufacturer
	Keep in a cool, dry place
	Caution, read instructions

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