Luxatemp-Glaze & Bond

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Luxatemp-Glaze & Bond is a light-curing, 1-component material based on multi-functional methacrylates for surface sealing and repair of temporary crowns and bridges, as well as for the pre-treatment of indirect plastic restorations.

Indications

- Surface sealer and repair material for temporary crowns and bridges
- As a primer/bonding agent for indirect restorations of PEEK, composite, and PMMA

Contraindications

- ▶ Not to be used if allergies are known to exist with components (especially methyl methacrylate).

Notes for use

- ▶ Light-curing units should have an output of 450 nm and should be checked regularly. Place the light as close as possible to the material. The light intensity should be a minimum of 400 mW/cm².
- ▶ Luxatemp-Glaze & Bond binds to teeth and films made of polycarbonate or poly acetate.
- ▶ Do not use Luxatemp-Glaze & Bond as a dentine or enamel adhesive.
- ▶ Avoid contact with materials containing eugenol.
- ▶ Close the bottle immediately after use.

Recommended use

Surface coating

1. Apply a thin layer of the material without entrapped air bubbles to the required surfaces.
2. Allow the material to act for 20 s and then cure each element for 20 s with a halogen or LED lamp (light intensity > 1000 mW/cm²). Alternatively cure the entire restoration for at least 90 s in the light furnace.

Note: To ensure that no oxygen inhibition layer is formed, the material should be light-cured within 2:00 min after being applied.

Adhesion for repair purposes

1. Roughen the surface to be repaired with a cutter.
2. Apply a thin layer of the material without entrapped air bubbles to the required surfaces.
3. Expose for 20 s with a halogen lamp or ≥10 s with an LED lamp (light intensity > 1000 mW/cm²).
4. Apply freshly mixed temporary material, e.g. Luxatemp Star, or a light-curing composite, e.g. LuxaFlow Star, to the relevant surfaces. Observe the manufacturer's instructions for use for the material used.

Bonding agent for restorations

1. Prepare the restoration according to manufacturer's specifications.
2. Apply a thin layer of the material without entrapped air bubbles to the required surfaces. Expose for 20 s with a halogen lamp or ≥10 s with an LED lamp (light intensity >1000mW/cm²).
3. Lute restoration with a suitable adhesive dual- or self-curing composite luting cement according to the manufacturer's instructions.

Storage and shelf life

- ▶ Store at 2–25 °C (36–77 °F) in a dry place.
- ▶ If stored in the refrigerator (2–8 °C / 36–46 °F), prior to applying bring the material to room temperature (15–25 °C / 59–77 °F).
- ▶ Do not use after the expiration date.

Composition

Multifunctional acrylates, methyl methacrylate, catalysts, stabilizers, additives

Packaging

Package	
1 bottle @ 5 ml liquid, 25 brushes	REF 211927
Accessories	
50 brushes	REF 110041

Keep out of the reach of children! For dental use only!

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