# Orthodontic Instruments

## Orthodontic Pliers and Orthodontic Instruments

**ErgoDenta ApS**

Hørskættø 14, 2630 Taastrup, Denmark.

Tel: +4555241100

Email: info@ergodenta.dk

Web: www.ergodenta.dk

en-English

### **Attention**

The operator/user must carefully read and understand this manual thoroughly to keep the product performance durable and reliable for defined life expectancy. After opening the packaging, first of all it is necessary to check the component(s) against the standard configuration. Check that these are all present and in perfect condition.

### **Intended Use:**

These instruments are Intended to assist in various orthodontic procedures, including the grasping, bending, cutting, adjusting, placement, and removal of orthodontic wires, brackets, and bands.

The Instrument Specific Intended Use along with their variants is given below:

### **Orthodontic Pliers:**

Orthodontic pliers are intended to be used for precise manipulation of orthodontic appliances, ensuring effective treatment outcomes.

### **Variants:**

**Nose-Pliers:** For bending and adjusting orthodontic wires.

**Universal Plier:** Versatile plier for various orthodontic applications like wire bending and adjustments.

**Aderer:** For delicate bending and forming of orthodontic wires.

**Weingart:** Helps in grasping and placing archwires securely in brackets.

**Schwarz:** For shaping and adjusting orthodontic wires.

**Adams:** Primarily for adjusting Adams clasps on removable appliances.

**Angle:** For precise wire bending and shaping.

**Tweed:** For torquing and bending orthodontic wires.

**Langback:** For forming and adjusting orthodontic wires.

**Ligature Cutter:** For cutting ligature wires and elastics in orthodontic treatments.

**Mathieu Ortho Plier:** For placing and securing elastics and ligatures.

**HOW:** Assists in placing and removing archwires in orthodontic brackets.

**Bracket Removing Plier:** For safe removal of orthodontic brackets.

**Adhesive Removal** – For scraping off residual adhesive after bracket removal.

**Cinch Back:** For bending distal ends of archwires to secure them.

**Buccal Tube Holder:** Helps in holding and placing buccal tubes on molars.

**Mathieu Orthodontic:** For handling ligatures and elastics in orthodontic procedures.

**Ribbon-tweed arch-short:** For bending and forming archwires in orthodontic treatments.

**Three beak:** Provides a secure grip with three distinct beak tips for intricate wire adjustments.

**Posterior band remover-long:** For removing posterior bands from molars.

**De la Rosa band contouring:** Helps in contouring and adjusting orthodontic bands.

**Separating module:** For placing elastic separators between teeth before band placement.

**Pin and ligature cutter:** For cutting ligatures and small-diameter orthodontic wires.

**Micro ligature cutter:** Designed for precision cutting of fine ligature wires.

**Distal end cutter:** Cuts and holds excess archwire behind the molars.

**Hard wire cutter:** For cutting thick orthodontic wires.

**Wire cutter plier:** Assists in cutting orthodontic wires and ligatures.

### **Orthodontic Instruments:**

These instruments assist in the placement, adjustment, and removal of orthodontic appliances.

### **Variants:**

**Band Pusher:** For pushing and adjusting orthodontic bands into their final position.

**Band Seater:** Assists in seating and positioning orthodontic bands securely onto teeth by applying controlled pressure.

**Boon's Gauge (With Lid):** for determining the ideal archwire size and dimensions, ensuring proper fit and alignment.

**Ortho Ligature Tucker:** Helps in tucking and securing ligatures around brackets.

**Ortho Archwire Tucker:** for securing archwires into bracket slots.

 Our Orthodontic Instruments and Pliers are intended to be used by qualified dental professionals only in a professional dental facility.

**Patient population:**

These instruments are intended for use on adult and child patients undergoing orthodontic treatment. The selection of the appropriate instrument type and size must be based on the patient's dental anatomy and the specific orthodontic procedure requirements.

**For Available Configurations/Models please visit the Orthodontics section in our Catalogue**

**Features**

**Orthodontic Pliers:**

- Manufactured using corrosion-resistant, medical-grade stainless steel for durability and longevity.
- Precision-aligned tips and ergonomic handles minimize user fatigue and enhance control.
- Textured or coated handles ensure a secure grip, even when wearing gloves.
- Designed to withstand repeated sterilization cycles without structural compromise.
- Finely machined and aligned jaws enable accurate bending, cutting, and crimping.

**Orthodontic Instruments:**

- Heat-treated, high-strength stainless steel ensures optimal hardness and wear resistance.
- Precisely crafted hinge mechanisms enable smooth and controlled movement.
- Electropolished surfaces enhance corrosion resistance and reduce microbial contamination.

**Operating instructions**

⚠ These instruments are non-sterile and must be cleaned and sterilized before first use.

⚠ Always check if the reprocessed instruments are in perfect condition and do not use the instruments if any kind of deterioration, damage or deformation is observed.

**Orthodontic Pliers:**

- Select the appropriate type and size of pliers based on the orthodontic procedure.
- Ensure the pliers' joint and working ends are properly aligned and functional.
- Hold the pliers with a firm but controlled grip for precision in handling wires, bands, or brackets.
- Apply gentle and steady pressure to avoid damaging wires, brackets, or patient soft tissues.

- Avoid excessive force, which can lead to distortion of orthodontic components or patient discomfort.

### **Orthodontic Instruments:**

- Handle with a secure grip to ensure precision during application.
- Use the instrument as per standard orthodontic techniques to adjust, place, or remove orthodontic appliances.

### **Adverse Event**

Any serious incident that has occurred concerning the device should be reported to the manufacturer and to the respective competent authority of the Member State in which the incident has occurred.

### **Cleaning & Sterilization**

#### **Cleaning**

⚠ Clean the device with care after each application to avoid the drying of contaminants.

#### **Initial Cleaning**

1. Remove any visible debris from the instrument using moisten a cotton swab or cloth with isopropyl or ethyl alcohol and gently wiping the surface of the product.
2. Use a soft brush and enzymatic cleaner to remove contaminants thoroughly.

#### **Manual Cleaning**

3. Prepare a solution of neutral pH enzymatic detergent according to the manufacturer's instructions.
4. Immerse in a pre-soak enzymatic cleaner solution
5. Thoroughly clean all surfaces of the Instrument using a soft brush or cloth. Pay close attention to any crevices or hard-to-reach areas.
6. Rinse thoroughly under lukewarm running tap water for a minimum of 30 seconds to remove all traces of detergent.

### **Automated Cleaning**

1. Always clean in an appropriate box or cassette and ensure that they are properly positioned within the unit.
2. Always follow the equipment manufacturer's instructions and ensure the instruments are compatible with the cleaning system.
3. Visually inspect for cleanliness, If any contamination is visible, repeat the cleaning steps until required cleanliness is achieved.
4. Dry with clean, lint free cloth or filtered compressed air until there is no visible moisture. Follow STEAM STERLIZATION PROCEDURE.

 Ultrasonic cleaning may be performed in accordance with hospital or facility protocols. Ensure proper rinsing and drying after the process to prevent residue buildup.

### **Sterilization**

Note: Ensure that the instruments are completely dry before sterilization.

### **Steam Autoclaving**

1. Place the instrument in an approved sterilization pouch or wrap.
2. Sterilize using a steam autoclave following a validated sterilization cycle according to ISO 17665-1 and your facility's protocols. A typical recommended cycle is 132-134°C (270-273°F) for 4 minutes.
3. Verify no residual moisture remains after sterilization to prevent corrosion and ensure longevity.
4. Ensure that pliers hinges are lubricated with an appropriate medical-grade lubricant after sterilization to maintain smooth operation.
5. Allow the instruments to cool and dry completely before handling. Store in a clean, dry, and controlled environment until ready for use.
6. Follow the pouch/wrap manufacturer instructions for storage conditions and maximum storage time.
7. For pliers, periodically check the hinge movement to ensure continued smooth function.

### **Important Notes**

- Always Inspect the instruments for any signs of deterioration, breakage and bent etc. after each sterilization cycle. Discard if damage is observed.
- Always follow your facility's infection control protocols.

- Refer to your autoclave manufacturer's instructions for proper operation and maintenance. The Sterilization equipment must be validated by the hospital and or sterilization equipment manufacturers.
- The above-mentioned cleaning and sterilization guidelines, provided by manufacturers, are intended as procedures compatible with specific materials. Sterilization must be performed according to the Hospital/Clinic approved protocol.
- The responsibility for ensuring proper sterilization of instruments lies with the user. The effectiveness of sterilization depends on validated cleaning, packaging, and sterilization procedures carried out at the facility.


### **Disposal**

Our Orthodontic Instruments and Pliers contain no hazardous materials. However, these must be cleaned/sterilized as per defined procedures before disposition. Please follow the local and national regulations or healthcare facilities' defined disposal and waste management policies.

### **Disclaimer**

The products must be used, reprocessed, and maintained strictly in accordance with the instructions provided above. Any deviation from these guidelines by the dental professional or user is undertaken at the user's sole risk. The Manufacturer will not accept any requests for refunds or exchanges under warranty for products that have not been handled and reprocessed in compliance with these instructions.

### **Explanation of utilized symbols**

	Medical Device Compliant with Regulation (EU) 2017/745
	Consult instructions for use
	Medical Device
	Product Code/Catalogue Number
	Model/Variant
	Lot Number
	The Device is supplied non-sterile.
	Sterilizable in a steam sterilizer (autoclave) at 132-134°C (270-273°F)

	Manufacturing Date
	Manufacturer
	Keep in a cool, dry place
	Caution, read instructions

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