# **PHILIPS**

Sonic electric toothbrush with app

## DiamondClean 9000

Connected brushing made easy

Built-in pressure sensor

Smart brush head recognition

4 modes, 3 intensities



HX9914/57



# Whiter, healthier teeth for life

Whiter teeth in just 1 day\*

Get the ultimate clean with DiamondClean 9000. Experience improved oral health with superior cleaning & whitening and real time feedback.

## DiamondClean users experience improved oral health

- Removes up to 10x more plaque\* for a deep clean
- ${\boldsymbol{\cdot}}$  Our unique technology brings you a powerful yet gentle clean

## Let the DiamondClean 9000 adapt to your needs

- · Personalized brushing experience
- · Lets you know if you're pressing too hard
- BrushSync automatically selects the best mode for you
- · Always know when to replace your brush heads

## Stay on track with a personalized progress report

· Start and keep up healthy habits

## Enjoy a stylish and user-friendly design

• Simple to charge and style to boot

## Highlights

## Removes up to 10x more plaque\*



Know you're getting the deepest possible clean with The C3 Premium Plaque Control brush head. The soft, flexible bristles are designed to curve around the contours of your teeth, giving you 4x more surface contact\*\* and up to 10x more plaque removal from hard-toreach spots

## **Philips Sonicare technology**



Powerful sonic vibrations whip up your toothpaste into plaque fighting bubbles, and drive them deep between your teeth and along your gum line. You'll get a month's worth of manual brushing in just 2 minutes.

## Built-in pressure sensor



You might not notice if you're brushing too hard, but your DiamondClean 9000 will. If you need to ease up, the toothbrush will make a pulsing sound. It's a heads up to let your brush head do the work. 7 out of 10 people found this feature helped them become a better brusher.

#### 4 modes, 3 intensities



The DiamondClean 9000 comes with Clean, White+, Gum Health and Deep Clean+ to take care of your brushing needs. Clean mode is for exceptional everyday cleaning, White+ is for removing stains, Gum Health provides a gentle yet effective clean for gums and Deep Clean+ gives you an invigorating deep clean. Three intensities allow you to choose between a higher setting to boost your clean and a lower one for more-sensitive mouths.

## BrushSync mode pairing



Smart brush heads ensure you're using the right mode and intensity for the best possible clean. For an example, say you're using the W3 Premium White brush head. Your DiamondClean 9000's BrushSync technology will automatically sync your brush head with the White+ mode to help whiten your teeth.

## **BrushSync replacement reminder**



All brush heads wear out over time, so you'll want to keep an eye on yours to make sure you're still getting a great clean. Our

BrushSync technology tracks how long you've been using your brush head for, and how hard you've been brushing. The BrushSync replacement reminder on your handle and a short beep gives you a heads up when it's time to replace it.

#### **Progress Report**



The DiamondClean 9000 gives you the guidance you need to improve and maintain healthy brushing habits in between dentist checkups. Built-in smart sensors let you know when you're using too much pressure, and by connecting your brushing experience to the Sonicare app, a personalized Progress Report helps you stay on track, to see how much you've improved over time.

## Long-life battery



Your DiamondClean 9000 will last for up to two weeks on a single charge with regular use.

## Specifications

#### Modes

3 intensities: High, Medium, Low Clean: For exceptional everyday clean Deep Clean+: For an invigorating deep clean Gum Health: Special attention to molars White+: To remove surface stains

#### **Items included**

Handle: 2 DiamondClean 9000 Brush heads: 2 C3 Premium Plaque Defense

Glass charger: 1 Charger base: 2

## Connectivity

**Bluetooth**® wireless technology: Connected brushing app

#### Compatibility

Android compatibility: Android phones, Bluetooth 4.0 enabled tablets iOS compatibility: iPad 3rd Gen or higher, iPhone 4S or higher, with iOS7 or higher, with iOS7 operational system

## **Design and finishing**

Color: Black, Rose gold

#### **Cleaning performance**

Performance: Removes up to 10x more plaque\*

Pressure feedback: Vibrates handle to alert

user

Timer: BrushPacer and SmarTimer

## **Smart sensor technology**

Pressure sensor: Alerts when brushing too hard

BrushSync Replacement Reminder: Always

know when to, replace brush heads

BrushSync Technology: smart brush head,

Connects smart handle and

#### Ease of use

Handle compatibility: Easy click-on brush

heads

Replacement reminder: To always ensure best

results, reminder icon lights up Handle: Slim ergonomic design

## **Technical specifications**

Battery: Rechargeable Battery type: Lithium ION

Operating time (full to empty): 14 days\*\*\*

## **Power**

Voltage: 110-220 V

## **Service**

Warranty: 2-year limited warranty

© 2021 Koninklijke Philips N.V. All Rights reserved.

Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. or their respective owners.

Issue date 2021-06-14 Version: 4.1.1



 $<sup>^{\</sup>ast}$   $^{\ast\ast}$  based on two periods of two-minute brushings per day