elements[™] IC

Obturation System

The new cordless elements[™]IC is the latest complete and reliable obturation system that combines a Downpack heat source with a Backfill extruder.

Get excellent, predictable results faster and easier: rely on its longer lasting battery (compared to elementsfree)*, inductive charging, and a comfortable, ergonomic design.

	Main competitor	elements™IC
Battery		
Battery Life	Up to 5 procedures	★ 16 canals (at devices default temperature, using recommended settings)
Battery Replacement	Battery not removable by the user	\star Battery removable with a screw driver
Battery Charge Time	5 hours	★ The battery fully charges in less than 3 hours
Connection		
Connection Base HPs	Pogo pins	★ Inductive Charging
Handpieces		
Gutta Percha Level Indicator	Yes 3 levels	★ Yes - 5 levels
Backfill length	21.5 cm	★ 20.5 cm
Grip	Notergonomic	📌 Ergonomic grip
Cartridge		
Cartridge Lock In	Difficult to remove after you bend it	🜟 Very easy to insert /remove

*Data on File

System and accessories order information and specifications

Reorder Code		
973-0600-TYPEF		
973-0600-TYPEG		
973-0610		
973-0612		
973-0615		
973-0616 -TypeF		
973-0616 -TypeG		
973-0620		

Description elementsIC Obturation System (EU) elementsIC Obturation System (UK) Dual Charger, elementsIC Single Charger, elementsIC Transformer elementsIC Power Cord, Type F (EU) Power Cord, Type G (UK) Battery, elementsIC

Dental Excellence in every area



Practice equipment

KaVo treatment units and lights, dental chairs, patient communication system, dental microscope and additional operatory accessories

Instruments

Dental straight handpieces and contra-angles, turbines, air polishing systems and small equipment for all application areas including diagnosis, prophylaxis, restorative, surgery, endodontics and instrument care



Imaging

Intraoral X-ray equipment, sensors and imaging plate systems, panoramic and cephalometric in combination with CBCT, as well as dedicated CBCT devices for every indication in dentistry

Endodontics

Restoratives

Solutions for diagnosis. shaping irrigation and obturation/filling.Plus suitable endodontic and ultrasonic products

Prevention



Prophy pastes, cups and brushes, dental lip retractor and trays, splash and X-ray protection, barriers for dental practice, X-ray holders, perio instruments, dental air/water syringe tips, sterile examination instruments, surface disinfectants, dental instrument reprocessing, dental office organization, and implant prophylaxis

The products, features and services shown and described in this catalogue are not available in all countries. All specifications were correct at the time of publication. KaVo Dental GmbH and KerrHawe SA assume no liability for deviations in colour or form from the illustrations, mistakes or printing errors and reserve the right to make changes to the brochures at any time. Reprinting, even just of excerpts, is only permitted with the consent of KaVo Dental GmbH and KerrHawe SA.

KaVo Dental GmbH | Bismarckring 39 | 88400 Biberach | Germany KerrHawe SA | Via Strecce 4 | 6934 Bioggio | Switzerland



Dental bonding agents, composites, curing lights, dental impression materials, permanent dental cements, liners and matrix systems as well as products for finishing and polishing as well as for dental restorations and alloys



Kavo Kerr transforms your Endo and takes it to the next level Just Endo it! Now you are In Control







Introducing the latest cordless elements™IC obturation system, the reliable companion for accurate, fast and excellent clinical outcomes.



Downpack

The Downpack device provides fast heating of the heat plugger with precisely controlled temperature and timing, making it suitable for singlemotion downpack obturation of the apical portion of the root canal.

- **Total control of the procedure:** adjustable heating range from 140°C to 400°C
- Ready when you are: heats to 200°C in less than 0.5 seconds
- Comfort without compromises: 360° activation ring

Downpack order information and specifications

Reorder Code 973-0602-TYPEF 973-0602-TYPEG Description Downpack Unit, elementsIC (EU) Downpack Unit, elementsIC (UK)

Downpack and Backfill common features

- Same screen
- Same digital interface
- Same controls
- Organic Light Emitting Diode display (OLED)

Backfill



- Same, easy to replace battery
- Temperature display
- Volume control: 6 levels
- Preset temperatures for gutta percha and custom setting