



Blueprint® Xcreme Dust-Free Alginate Impression Material

- Smooth compact surface
- Excellent elastic recovery
- 5 days of dimensional stability

Blueprint® Xcreme

Dust-Free Alginate Impression Material

For many years, Blueprint®-Alginates have been the backbone in the clinical procedures of countless dental professionals all over the world. Blueprint® Xcreme is a state of the art Alginate for the needs of today's dental professionals. It offers ...

- **Smooth compact surface through extra fast water absorption**
- **Excellent elastic recovery (98%)**
- **5 day dimensional stability with minimal water loss**



Clinical Benefit

Blueprint® Xcreme impressions provide excellent quality which is achieved easily and conveniently. Its dimensional stability contributes to getting a high-quality stone model out of a high quality impression.

Order Information

Blueprint® Xcreme Refill Pack 605.01.100
 2 x 500 g Blueprint® Xcreme
 1 x Storage Tub
 1 x Spoon & Cup

Blueprint® Xcreme Eco Pack 605.01.108
 12 x 500 g Blueprint® Xcreme
 1 x Storage Tub
 1 x Spoon & Cup

Blueprint® Xcreme Measures 605.01.998
 1 x Spoon & Cup



Fix Adhesive
 Tray Adhesive for Alginate
 Impression Materials
 1 x 200 g Aerosol
 1 x 30 ml Liquid

616.01.001
616.01.004



Fix Solvent
 Cleanser for Removal
 of Fix Adhesives
 1 x 200 g Aerosol

616.01.101



For more information visit www.dentsplysirona.com.

Find us on  www.facebook.com/dentsplysirona.restorative

K79200040-03
 © Dentsply Sirona 09/2018