

let's keep connected



Polydentia SA
via Cantonale 47
6805 Mezzovico-Vira
Switzerland

phone: +41 91 946 29 48
fax: +41 91 946 32 03

info@polydentia.ch
www.polydentia.ch



Vista-Tec ultra light

protective
face- and
eyeshields



REV 03 09 2017



polydentia
swiss manufacture

Vista-Tec ultra light

LATEX
FREE

reusable
shields

frames
available in
six translucent
colors



autoclavable
frame



Vista-Tec Faceshield

Feather-light
protective faceshield



Vista-Tec provides full facial protection and an excellent optical view. It is comfortable to wear and elegantly designed; in addition, it does not hinder communication with the patient.

The protective shields, which are made of the highest quality plastic, are scratch resistant, antiglare and antifog. They can be used several times as they can be cleaned easily with mild soap and water. Shields are also available separately as refills.

The frame is lightweight, flexible and can be worn comfortably all day. It is available in six translucent colors. Available as well in the black autoclavable and thermo-disinfectable version.

The clip integrated into the frame allows easy fixing of the protective shield. Once inserted on the frame, the shield is angle and height adjustable.

Vista-Tec ultra light Eco: 1 frame, 5 shields

white	REF 5621/W
yellow	REF 5621/Y
blue	REF 5621/B
transparent	REF 5621/T
pink	REF 5621/P
green	REF 5621/G

Autoclavable Vista-Tec: 1 frame


black	REF 5613
-------	----------

Vista-Tec ultra light Classic: 1 frame, 10 shields

white	REF 5624/W
yellow	REF 5624/Y
blue	REF 5624/B
transparent	REF 5624/T
pink	REF 5624/P
green	REF 5624/G

Refills:

5 Faceshields	REF 5617
10 Faceshields	REF 5616



Clipon & L-Clipon

Protective shields for
glasses and loupes



Clipon - The plastic clips to fit Vista-Tec shields on prescription spectacles, supplied with classical faceshields. For L-Clipon, the plastic clips are coming with L-Faceshields, a larger version of the shields fitting over loupes.

Clipon, for spectacle wearers

2 Clipon, 5 shields	REF 5622
5 faceshields	REF 5617
10 faceshields	REF 5616

L-Clipon, for loupes

2 Clipon, 5 L-shields for loupes	REF 5630
10 L-faceshields	REF 5632



Vista-Tec Eyeshield

Feather-light
protective eyeshield



The shields are also available in the eyewear version suitable for patients as well. The special design of these shields provides comfortable and complete protection when worn with surgical masks. The lightweight frame is flexible and comfortable thus when worn, no pressure is exerted on the head, ears and bridge of the nose. It is available in white and five delicate translucent colors: transparent, yellow, blue, pink and green.

Vista-Tec ultra light Eco: 1 frame, 5 eyeshields

white	REF 5628/W
yellow	REF 5628/Y
blue	REF 5628/B
transparent	REF 5628/T
pink	REF 5628/P
green	REF 5628/G

Refills:

5 Eyeshields	REF 5618
10 Eyeshields	REF 5619

Autoclavable Vista-Tec: 1 frame

black	REF 5613
-------	----------



Vista-Tec Orange shields

Perfect protection for
high-intensity light



Protecting your eyes from dangerous blue light frequencies coming from curing lights. The high quality of these orange-tinted shields allows an excellent optical view, while keeping the eyes safe. These shields are easily fixed on Vista-Tec frames or Clipon to assure a comfortable, lightweight and effective barrier against high-intensity light. Available as full-face shield and as eyeshield, which is also suitable to protect patients.

Vista-Tec Orange shields

3 Orange faceshields	REF 5617/O
3 Orange eyeshields	REF 5618/O