

# **SOPIRA® Carpule needles**The right needle for any indication

Experience the ease of a gentle injection - your benefits at a glance:

- » Extensive range of dental needles for various dental anaesthesia techniques
- » Less initial pain due to the triple-bevel cut
- » Silicone coating ensures optimal penetration of tissue
- » Safe handling thanks to high elasticity
- » A variety of models to accommodate your surgery's every need: exceptionally thin-walled (Free Flow), safely packed (blistered) or in an extra short cut, especially for intraligamentary anaesthesia.
- » Fast and sure identification of the various needles with the new rainbow colour coding system
- » Compatible with all common dental syringes.

The right needle for any indication.





Giving a hand to oral health.



### **SOPIRA®**

The right needle for any indication.

## SOPIRA® Carpule needles: Gentle injection thanks to workmanship of the most outstanding quality.

The best injection is the one your patient doesn't feel. Our tissue-friendly needles are characterised by first-class steel and an innovative cut. They offer both dentist and patient a pleasant and comfortable treatment under ideal conditions.



#### SOPIRA Carpule needles

- Less initial pain due to three-faceted cut.
- Silicone coating ensures optimal gliding through the tissue.
- High elasticity ensures stable handling even with very long needles.
- The new colour-coding system for quick and reliable identification of the various needles.

Intraseptal
27G 08 mm



#### SOPIRA Carpule Free Flow needles

- Interior diameter as much as 50 % larger thanks to extra-thin walls.
- Low injection resistance for smoother handling.
- More comfortable for patients because the injection solution penetrates in a way that is even gentler on tissue.

Infiltrative	Intraseptal	Mandibular
27G 25 mm	_	27G 38 mm
30G 25 mm		

#### SOPIRA Carpule needles, blistered

• Individually sealed needles for even better safety, especially during surgery.



Infilt	rative
27G	25 mm
30G	$16\mathrm{mm}$
30G	21mm
30G	25 mm

Intraseptal 27G 08 mm Mandibular 27G 38 mm

Mandibular

27G 38 mm





Left picture: SOPIRA Carpule needles (from left to right): long bevel, Free Flow, short bevel Right picture: Free Flow needle with larger interior diameter