


UNIVERSALLY
STRONG.
UNIVERSALLY
SIMPLE.





SCAN TO
LEARN MORE!



CASE STUDY

Courtesy of Adamo Notarantonio, DDS, AACD, FICOI



ORDERING INFORMATION

DUO-LINK UNIVERSAL SYSTEM KIT..... A-19620K

1 Dual-Syringe DUO-LINK UNIVERSAL Shade: UNIVERSAL (8g), 1 Bottle ALL-BOND UNIVERSAL (6ml),
1 Bottle Z-PRIME PLUS (2ml), 1 Bottle Porcelain Primer (3ml), Accessories, Instructions

1 DUAL-SYRINGE PACKAGE

1 Dual-syringe (8g), 15 Mixing Tips and Intraoral Tips, Instructions

UNIVERSAL A-19030P

MILKY WHITE A-197MWP

ALL-BOND UNIVERSAL (6ml Bottle) B-7202P

Z-PRIME PLUS (2ml Bottle) B-6002P

ZIRCLEAR® (5g Syringe) B-7351P

PORCELAIN PRIMER (10ml Bottle) B-2221P

Call us! We're here to help:

1-800-247-3368

www.bisco.com



DUO-LINK UNIVERSAL™



Rx Only

DUO-LINK UNIVERSAL™

The adhesive resin cement is specially formulated for cementation of ALL† indirect restorations. It is available in two esthetic shades, Universal and Milky White. Duo-Link Universal is intended for use with adhesives designed for compatibility with all dental materials, including all BISCO adhesives.



THE PERFECT FIT FOR YOUR OFFICE



Extremely easy clean-up
Formulated to allow for quick and easy removal of excess cement



Highly radiopaque
Easily identified on radiographs for quick and effective diagnosis



CAD/CAM restorations
Ideal for all chairside and lab-fabricated restorations



Universal for all cementation procedures†
Indicated for all crowns, bridges (fixed prosthesis), inlays, onlays, and posts/dowels, fabricated from metal, composite, porcelain, glass ceramics, zirconia, and alumina.



Easy to use
Auto-mix, dual-syringe provides a consistent mix for immediate delivery



High degree of conversion
High degree of conversion in both light- and self-cured modes ensures a strong, long lasting restoration



Two shades
Universal shade is a translucent A2 and the Milky White shade is more opaque and can be used to mask darker tooth preps

† It is recommended to use BISCO's CHOICE™ 2 for veneer cementation.

WHAT'S IN THE KIT?

The Duo-Link Universal kit is an adhesive cementation system. The combination of award-winning, innovative products provides high retention and durability of indirect substrates in a simple, universal, and easy-to-use kit.

DUO-LINK UNIVERSAL

This adhesive resin cement specially formulated for cementation of ALL† indirect restorations

ALL-BOND UNIVERSAL

Light cure, MDP containing, truly universal, single bottle adhesive

PORCELAIN PRIMER

Silane coupling agent for all glass ceramics surfaces


Z-PRIME™ PLUS

Single component and MDP containing primer for zirconia and metal surfaces




SEE WHAT THE DATA SAYS


SBS Zirconia (Sandblasted)



SBS e.max (Etched)



SBS Dentin (Total Etch)




*All cements were used in light-cured mode. Trademarks are property of their respective owners.

DUO-LINK UNIVERSAL PREP AND CEMENTATION STEP BY STEP


1

RESTORATION TREATMENT**



2

ETCHING



3

ADHESIVE APPLICATION



4

CEMENTATION



** Refer to complete instructions for use for detailed technique and instructions.