

Tubulitec Primer and Tubulitec Liner[®]

- A COMPLETE SYSTEM FOR CAVITY INSULATION

Tubulitec Primer and Tubulitec Liner are a complete system for cavity insulation under all kinds of filling materials.

By first applying the dentin varnish, Tubulitec Primer, all dentinal tubular apertures are sealed in order to avoid micro leakage and effectively eliminate any microorganisms. The Tubulitec Primer will successfully adhere to the liquid Tubulitec Liner which is then applied as a second step to the interior of the prepared cavity before insertion of any kind of filling material.

In addition, Tubulitec Primer can be used as an antibacterial insulation under temporary crowns and bridges. With its multipurpose use Tubulitec Liner can be applied to protect the pulp of a tooth from chemical irritation.

Facts and benefits

- **Antibacterial primer** - Tubulitec Primer will kill all bacteria while sealing the dentinal tubular apertures
- **Universal varnish** - Tubulitec Primer can be used under all kinds of filling materials. Applied in seconds with a cotton pellet
- **Liquid cavity liner** - Tubulitec Liner, unlike pastes, can be applied to the walls of cavity preparation
- **No paste-like features** - As the composite begins to polymerize, Tubulitec Liner will not lift of where applied



Did you know:

Can be used under all kinds of filling materials. Applied in seconds with a cotton pellet eliminating any kind of bacteria infiltration.



Dental Therapeutics

phone: +46(0)8 716 81 41

email: info@therapeutics.se

web: www.therapeutics.se