


HARTMETALLFRÄSER | TUNGSTEN CARBIDE CUTTERS | FRESAS DE CARBURO PARA LABORATORIO

S


Standardverzahnt (Ref.-Nr. 175), blauer Ring, für wirkungsvolles Ausarbeiten aller harten Werkstoffe mit besonders glatten Oberflächen

Plain cut, medium (Ref.-No. 175), blue ring, for effective trimming of hard materials and for obtaining especially smooth surfaces

Dentado simple, medio (Ref.-No. 175), anillo azul, para desbaste eficaz sobre todos los materiales duros con superficies especialmente lisas



| Schaft / Shank | Ref.-Nr. / Ref.-No. | HM 21 | HM 23LR | HM 71 | HM 71 | HM 71 | HM 71 | HM 71 | HM 71 | HM 72 | HM 73 | HM 75 | HM 75 |
|----------------|---------------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|-------|
| | | 060 | 023 | 025 | 027 | 031 | 040 | 050 | 060 | 060 | 060 | 070 | |
| | L mm | 13,7 | 11,5 | 0,0 | 0,0 | 0,0 | 0,0 | 0,0 | 13,7 | 10,0 | 12,0 | 13,2 | |
| HP | 500 104 ... | ... 107 175 | ... 199 175 | ... 001 175 | ... 001 175 | ... 001 175 | ... 001 175 | ... 001 175 | ... 137 175 | ... 277 175 | ... 260 175 | ... 260 175 | |
| | | 2 | 2 | 2 | 2 | 2 | 2 | 2 | 2 | 2 | 2 | 2 | |
| RA | 500 204 ... | | | | | | | | | | | ... 260 175 | |
| | | | | | | | | | | | | 2 | |



| Schaft / Shank | Ref.-Nr. / Ref.-No. | HM 77 | HM 77 | HM 78 | HM 79 | HM 79 | HM 79 | HM 251 | HM 251 | HM 416 | HM 489 |
|----------------|---------------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|
| | | 050 | 060 | 060 | 040 | 045 | 070 | 060 | 070 | 060 | 023 |
| | L mm | 10,7 | 11,5 | 12,7 | 14,2 | 12,7 | 14,7 | 14,5 | 14,7 | 14,7 | 15,0 |
| HP | 500 104 ... | ... 237 175 | ... 237 175 | ... 257 175 | ... 194 175 | ... 194 175 | ... 194 175 | ... 274 175 | ... 274 175 | ... 201 175 | ... 187 175 |
| | | 2 | 2 | 2 | 2 | 2 | 2 | 2 | 2 | 2 | 2 |
| RA | 500 204 ... | | | | | ... 194 175 | | | | | |
| | | | | | | 2 | | | | | |




G

Grobverzahnt (Ref.-Nr. 215), grüner Ring, für schnelles, effektives Arbeiten auf Kunststoff und Gips

Plain cut, coarse (Ref.-No. 215), green ring, for rapid and effective trimming of acrylics and plaster

Dentado simple, grueso (Ref.-No. 215), anillo verde, para desbaste rápido y eficaz sobre acrílico y yeso



| Schaft / Shank | Ref.-Nr. / Ref.-No. | HM 75G | HM 75G | HM 79G | HM 79G | HM 251G |
|----------------|---------------------|-------------|-------------|-------------|-------------|-------------|
| | | 060 | 070 | 045 | 060 | 060 |
| | L mm | 12,0 | 13,2 | 12,7 | 12,7 | 14,7 |
| HP | 500 104 ... | ... 260 215 | ... 260 215 | ... 194 215 | ... 194 215 | ... 274 215 |
| | | 2 | 2 | 2 | 2 | 2 |

