

TECHNICAL DATA SHEET



Description

Reference

50004 – SV50A:
Safeview Eyewear Assembled Glasses/Assorted Colors
50002 – SV100F:
Safeview Eyewear Frames/Assorted Colors
50001 – SV250L:
Safeview Eyewear Replacement Lenses

Eyewear protects the eyes from splashing, spraying or splattering of liquids that may be irritating or injurious to the eyes or potentially infectious to the user.

Lightweight for added comfort in use.

Easily worn over most prescription glasses.

Main Materials

Lenses film: Clear Polyester

Frames: Polypropylene

Does not contain Natural rubber latex.

Does not contain DEHP.

Sterilisation

Products are non-sterile.

Packaging

50004: Shipping case of 50 units.

10 assemble glasses are placed within 1 sealed polyethylene bag, 5 sealed polyethylene bags are placed within 1 shipping case.

Bar coding: GS1-128 symbology, linear, on shipping case.

50002: Shipping case of 100 units.

10 reusable frames within 1 sealed polyethylene bag,

10 sealed polyethylene bags are placed within 1 shipping case.

Bar coding: GS1-128 symbology, linear, on shipping case.

50001: Shipping case of 250 units.

25 lenses are placed within 1 envelope, 10 envelopes are placed within 1 shipping case,

Bar coding: GS1-128 symbology, linear, on shipping case.



Manufacturing

Products are manufactured in the USA and Mexico.

Regulatory Information

Products CE marked as per 89/686/EEC Directive on Personal Protective Equipment.

Class of the devices: PPE category 1

Storage Information

Store in a dry and cool place, away from intense sources of heat.

Keep as much as practicably possible in its shipper box.

Shelf Life

None Specified.