

X-Smart® Plus

Everything you like about X-Smart® with a plus




Everything you like about X-Smart® with a plus



Reciprocation

Single file root canal shaping with WaveOne Gold the latest reciprocating technology to shape the canal with confidence.





Continuous motion

The next generation of the endo gold standard: ProTaper Next*, the performance of ProTaper* Universal with added versatility and a shorter clinical sequence.







Technical features

- Miniature contra-angle head that can be adjusted in 6 positions
- On/Off button on the motor handpiece
- Reciprocating motion and continuous rotation
- Includes a file library with preset programs for WaveOne*, WaveOne* Gold, ProTaper Next*, ProTaper* Universal, ProGlider*, PathFile*, Gates, Reciproc* - 6 free programs for individual settings in continuous rotation
- Auto reverse rotation at preset torque limit
- Torque range: 0.6 4.0 Ncm
- Speed range: 250 1'200 rpm
- Operates on a rechargeable battery and when connected to the mains
 - Battery autonomy: approx. 2 hours continuous use
 - Full recharge in approx. 5 hours







X-Smart® Plus WaveOne® Gold



- X-Smart Plus Endo motor
- 5 blisters WaveOne Gold primary, 25mm (15 files)
- 4 blister WaveOne Gold assortment, 25mm (16 files)
- 1 blister ProGlider, 25mm (3 files)
- WaveOne Gold clinical procedure card

X-Smart[®] Plus ProTaper Next[®]



- X-Smart Plus Endo motor
- 10 blisters ProTaper Next assortment X1-X3, 25mm (30 files)
- 1 blister ProGlider, 25mm (3 files)
- ProTaper Next clinical procedure card

