

Philips Sonicare
DiamondClean 9000 Series
Power Toothbrush Special
Edition

Connected app tracks progress

Built-in pressure sensor
Smart brush head recognition
4 modes, 3 intensities

HX9911/84



An effortlessly white smile

Whiter teeth in just 3 days*

Start every day in style with advanced Sonicare technology for whiter teeth and oral health. Gentle on gums with a connected app for coaching and tracking. Choose from 4 modes and 3 intensities for a personalized brushing experience.

Whiter teeth & improved oral health

- Hard on plaque, gentle on your gums

Adapts to your needs

- 4 modes & 3 intensities for personalized brushing experience

Adapts to your needs

- Pressure sensor lets you know if you're pressing too hard
- BrushSync automatically selects the best mode for you

Stay on track with a personalized progress report

- Sonicare app for coaching and progress reports

Stylish, user-friendly design

- Stylish design for a premium brushing experience

PHILIPS
sonicare

Highlights

Whiter teeth in just 3 days*



Boost your smile's brightness with advanced Sonicare technology. Polishing bristles remove up to 100% more stain than a manual toothbrush for visibly whiter teeth in just 3 days*.

Gentle on your gums



The C3 Premium Plaque Control brush head is designed to give you your deepest clean yet. Soft, flexible sides and bristles perfectly curve around the contours of your teeth, giving you 4x more surface contact and helping target hard-to-reach spots***

Premium brushing experience



Our DiamondClean 9000 comes with a standout ombré design for a premium, ergonomic brushing experience that suits your style. Slim, compact travel case in a complementary color goes anywhere life takes you.

Advanced Sonicare technology



Powerful sonic vibrations move 62,000 times every minute to drive plaque-fighting bubbles between your teeth and along your gum line. Removes up to 10x more plaque than a manual toothbrush.

4 modes / 3 intensities



Choose from 4 brushing modes: Clean for exceptional everyday cleaning, White+ to brighten your smile, Gum Health for a gentle yet effective clean, and Deep Clean+ for an invigorating deep clean. Three intensity settings from Low to High let you boost performance or care for sensitive gums.

Built-in pressure sensor



You might not notice if you're brushing too hard, but your DiamondClean 9000 will. If you need to ease up, the toothbrush will make a pulsing sound. It's a heads up to let your brush head do the work.

BrushSync automatic pairing



Smart brush heads ensure you're using the right mode and intensity for the best possible clean. When you are using the C3 Premium Plaque Control brush head, your DiamondClean 9000's BrushSync technology will automatically sync your brush head with the Clean mode to help whiten your teeth.

Sonicare app



The DiamondClean 9000 gives you the guidance you need to improve and maintain healthy brushing habits in between dentist checkups. Built-in smart sensors let you know when you're using too much pressure, and by connecting your brushing experience to the Sonicare app, a personalized Progress Report helps you stay on track, to see how much you've improved over time.

Specifications

Design and finishing

- Color: Silk Pink to White Gradient

Connectivity

- Bluetooth® wireless technology: Connected brushing app

Power

- Voltage: 110-220 V

Technical specifications

- Battery: Rechargeable
- Operating time (full to empty): 14 days***
- Battery type: Lithium ION

Service

- Warranty: 2-year limited warranty

Compatibility

- Android compatibility: Android phones, Bluetooth 4.0 enabled tablets
- iOS compatibility: iPad 3rd Gen or higher, iPhone 4S or higher, with iOS7 or higher, with iOS7 operational system

Ease of use

- Handle compatibility: Easy click-on brush heads

- Handle: Slim ergonomic design
- Replacement reminder: To always ensure best results, reminder icon lights up

Items included

- Brush heads: 1 C3 Premium Plaque Control
- Handle: 1 DiamondClean 9000
- Travel case: 1
- Charger base and puck: 1

Cleaning performance

- Performance: Removes up to 10x more plaque*
- Timer: BrushPacer and SmarTimer
- Pressure feedback: Vibrates handle to alert user

Modes

- 3 intensities: High, Medium, Low
- Clean: For exceptional everyday clean
- Deep Clean+: For an invigorating deep clean
- Gum Health: Special attention to molars
- White+: To remove surface stains

Smart sensor technology

- BrushSync Replacement Reminder: Always know when to, replace brush heads
- Pressure sensor: Alerts when brushing too hard



Issue date 2022-05-25

Version: 11.11.0

© 2022 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com

* in White+ Mode using a leading whitening toothpaste
** Vs a manual toothbrush, actual results may vary
*** as compared to DiamondClean brush head