

# ***PULPDENT Corporation***



## **Pic-n-Stic**

**It's the extra fingers you've  
been waiting for**


# Pic-n-Stic

The original “handle for small objects” has numerous applications from dentistry to model making and placing hearing aid batteries.



Pic-n-Stic is 2” long 2 mm in diameter with an adhesive tip on one end.

**Pic-n-Stic is great for holding inlays and onlays while etching, placing silane, bonding agents and cement.**



***PULPDENT***

**Pic-n-Stic is perfect for holding and placing veneers.**



***PULPDENT***



# Pic-n-Stic

“The original handle for small objects”

Order Code: PIC Box of 60