

Enhance® Finishing System

Coarse when you need it, fine when you want it.

The Enhance Finishing System accomplishes the entire process of excess removal, contouring, and finishing simply by varying the pressure. Push harder to remove excess composite, ease it off for final results. In just one step the Enhance Finishing System delivers a smooth, contoured surface with a natural looking finish. It's all about pressure. **Finish in one step.**




Enhance® mini

Improved access and versatility

Fewer Steps. Less Time. Impressive Results.


More than 70% of dentists say that simplicity is one of the most important factors for a finishing and polishing system.¹ That's why we designed the Enhance Finishing System. Achieve the entire finishing process in one step. It doesn't get simpler than that!

Apply more pressure for excess removal and initial contouring; back it off for a naturally smooth surface. With no need to change instruments, you can stay focused on the procedure to achieve a better finish, faster.



Maximum gloss per system

After achieving a natural looking finish with Enhance, create a high-end luster with the Enhance® PoGo® diamond polishing system.



Ordering Information

Enhance® Finishing System Refills

40 Finishing Discs	624045
40 Finishing Cups	624055
40 Finishing Points	624065

Enhance®/Enhance® PoGo® Complete Kit

Enhance® Finisher	10 each: Discs, Cups, Points	662020
Enhance® PoGo® Polisher	10 each: Discs, Cups, Points	

Enhance® mini Intro Kit

Enhance® mini	20 each: Discs, Cups, Points	624385
---------------	------------------------------	---------------

Enhance® mini Finishing System Refills

40 Finishing Discs	624345
40 Finishing Cups	624355
40 Finishing Points	624365



For more information visit www.dentsplysirona.com

Find us on www.facebook.com/dentsplysirona.restorative



1. Exevia: Dental Shuttle Q4, 2016, Germany
 2. da Costa J, Ferracane JL (OHSU, Portland, OR)
 3. Not a registered trademark of Dentsply Sirona, Inc.
K79200260-00

© Dentsply Sirona Inc. 06/2017