

# Sendoline Endo Motor - maximum freedom of choice!

Sendoline Endo Motor is a cordless rechargeable motor for endodontics. The latest technology unites the simplicity of preset programs with the freedom to choose your own settings – maximum freedom of choice.



#### Access

The neat contra-angle is rotatable 360° allowing easier access, while the narrow neck and small head of the contra-angle do not obscure your view.

#### Choice

Sendoline Endo Motor can be used for both rotary and reciprocating systems, giving you maximum freedom of choice.

## Sendoline® Endo Motor



### With or without an apex locator – You choose! Sendoline Endo Motor can be used in 3 different ways:

- 1. Only as an endo motor
- 2. Endo motor with apex locator
- 3. Only as an apex locator.





Scan the QR code for more information about Sendoline Endo Motor

