

OptiBond FL®. Proven strength. New Unidose® delivery.



OptiBond/OptiBond FL wins REALITY'S 20th Anniversary Legacy Award, emblematic of extraordinary long-term clinical performance.



TOTAL-ETCH ADHESIVE

The best gets even better.

- **Versatile.** Two-step system for all direct restorations or indirect core buildup.
- **Easy-to-use.** One coat primer. One coat adhesive. Wet or dry prep.
- **Unique structural bond.** 48% fill load delivers superior bond strength.
- **Delivery options.** Choice of bottle or Unidose.
- **Proven long-term performance.** The gold standard in adhesive dentistry after more than a decade on the market.

Order Information

26684 OptiBond FL Syringe Kit – Contains 1 Bottle Prime (8 ml), 1 Bottle Adhesive (8 ml), 1 Syringe Kerr Gel Etchant (3 g), 10 Gel Etchant Syringe Tips (disposable), 25 Mixing Wells (disposable), 50 Kerr Applicators (disposable) and Technique Card

33352 OptiBond FL Unidose Kit – Contains 50 Unidose Packs (15.5 g total) – Each pack contains 1 Unidose Prime (yellow), 1 Unidose Adhesive (black) and 1 Applicator

25881 FL Prime Refill (Bottle #1 – 8 ml)

25882 FL Adhesive Refill (Bottle #2 – 8 ml)

Explore the science of Kerr at www.KerrDental.com • 800.KERR.123

