# SELECT HV. ETCH

## 35% Phosphoric Acid Etchant

SELECT HV ETCH is a 35% high viscosity phosphoric etch. SELECT HV ETCH is used to condition the tooth structure before bonding adhesives, composites, or sealants. SELECT HV ETCH is specially formulated for maximum working and handling, pin-point placement performance, and eliminating run-on onto the occlusal dentin surface.

- High viscosity, ideal for enamel etching
- Thixotropic
- Blue in color for easy visualization and contrast
- Washes off easily without leaving residue
- Bulk syringe dispensing



SELECT HV ETCH delivers smooth, non-stringy material enabling it to perform well during the "hybrid technique" – etching enamel margins without etching dentin. This is desirable when etching enamel in conjunction with applying a self-etch adhesive, or for the immediate dentin sealing procedure.



SELECT HV ETCH indicating pin point placement on enamel.

Dentistry courtesy of Dr. Michael Morgan

**To Order Call:** 1-847-534-6000











MC-3258SEI Rx Only





## 35% Phosphoric Acid Etchant



### **INDICATIONS FOR USE:**

- 1. Etching dentin and enamel
  - Time: 15 seconds
- 2. Selective enamel etching
  - Time: 15 seconds
- 3. Cleaning agent on dental restorative materials
  - Time: 30 seconds

#### ORDERING INFORMATION

#### KIT CONTENTS

SELECT HV ETCH

1 Syringe SELECT HV ETCH (30ml), 30 Disposable Syringes,

30 Disposable Tips, Instructions

#### **REFILLS**

#### SELECT HV ETCH

1 Bulk Syringe (30ml), Instructions E-59060P
4 Syringes (5g ea.), 50 Disposable Syringe Tips, Instructions E-59010P

EMPTY SYRINGE ACCESSORY PACK

30 Disposable Syringes, 30 Disposable Syringe Tips, Instructions X-80580P

50 Dark Blue Disposable Syringe Tips (22 Gauge) X-80608N